

# Hip Retractors

Symmetry Medical's Hip Retractors are designed for use in various surgical exposures to the hip, including minimally invasive approaches.

Typical Hohmann-type blade ends are combined with unique handle shapes for improved visualization and easy wound access. Each retractor thickness is optimized for ease of placement, visualization and stability.



## BENEFITS

- Typical Hohmann-type blade ends.
- Retractor bend angles remove handles from wound site for improved visualization.
- Optimized thickness helps wound integrity and stability.
- Can be used for various minimally invasive hip exposures.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Hip Retractor System	30-034-000
Hip Retractor, 25mm Blade	13-023-000
Hip Retractor, Blunt Tip	13-024-000
Offset Retractor, 32mm Blade	13-025-000
Long Hohmann	13-026-000
Long Hohmann, Curved, Short Blade	13-027-000
Long Hohmann, Curved, Long Blade	13-028-000
Long Hohmann, Narrow	13-029-000
Small Hohmann	13-030-000
Femoral Elevator	13-031-000
Hybrid Case	25-028-000



**SYMMETRY Medical**<sup>®</sup> Inc.

486 West 350 North, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.7306  
www.symmetrymedical.com