

Symmetry Medical Instrument Cassettes



The multifunctional Instrument Cassettes accommodate a wide variety of delicate instruments.

Instrument Cassettes are ideally suited for use with ophthalmic, ear/nose/throat, general dentistry, cosmetic and microsurgery instruments

BENEFITS

- Hinged lid provides easy use.
- Silicone mats gently secure instruments.
- Raised silicone (latex free) mats provide ventilation to facilitate drying.
- Silicone bars accommodate instrument diameters 2.5 - 8mm.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Shallow Instrument Cassette with Silicone Mat, 8 x 4.5 x .61 inches	7-6050
Deep Instrument Cassette with Silicone Mat, 8 x 4.5 x 1.2 inches	7-1250
6-Position Slotted Instrument Cassette with Silicone Bars, Holds up to 6 Hand Held Instruments 2.5 - 8mm Diameter, 8 x 4.5 x .61 inches	7-6000
12-Position Slotted Instrument Cassette with Silicone Bars, Holds up to 12 Hand Held Instruments 2.5 - 8mm Diameter, 8 x 4.5 x 1.2 inches	7-1200



SYMMETRY Medical^{INC}

253 Abby Road, Manchester, NH 03103, USA
T: +1.603.647.7822 Fax: +1.603.647.7877
www.symmetrymedical.com

ZI de la Pilaterie, 14, Rue de la Ladrie, BP 165, 59653 Villeneuve d'Ascq Cedex, France
T: +33.3.28.09.94.54 Fax: +33.3.20.72.31.04