

42.5 Degree Driver

The Symmetry Medical 42.5 Degree Driver* is an extension for drilling and reaming. The driver allows the surgeon access to difficult anatomical areas, specifically glenoid and acetabular applications.

BENEFITS

- 1/2" Driver is used in hip surgery to drill and tap angled holes.
- 3/4" Driver is used to prepare the glenoid surface for component seating.
- Available in multiple driver ends with various tool connections.
- The unique universal joint has a smooth 360 degree swivel end.

**The 1/2" Driver cannot be used in shoulder applications.*

Contact your local Sales Manager for more information.



42.5 Degree Driver

DESCRIPTION	PART #
1/2" Hudson End - Bayonet	09-078-000
1/2" Hudson End - Threaded	09-079-000
1/2" Tri-Shank End - Bayonet	09-080-000
1/2" Tri-Shank End - Threaded	09-081-000
1/2" A/O End - Bayonet	09-082-000
1/2" A/O End - Threaded	09-083-000
5/8" Hudson End - Bayonet	09-133-000
5/8" Hudson End - Threaded	09-134-000
5/8" Zimmer End - Bayonet	09-135-000
5/8" Zimmer End - Threaded	09-136-000
5/8" A/O End - Bayonet	09-137-000
5/8" A/O End - Threaded	09-138-000
3/4" Hudson End - Bayonet	09-084-000
3/4" Hudson End - Threaded	09-085-000
3/4" Tri-Shank End - Bayonet	09-086-000
3/4" Tri-Shank End - Threaded	09-087-000
3/4" A/O End - Bayonet	09-088-000
3/4" A/O End - Threaded	09-089-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com