

# Modular Calcar Planer

The Symmetry Medical Modular Calcar Planer provides a safe, effective and accurate cutting instrument for the removal of hard calcar bone.

The patented design provides easy disassembly and blade change in a matter of seconds. Surgeons will appreciate a constantly sharp blade producing repeatable, accurate cuts.



## BENEFITS

- Disposable blades provide sharp cuts for every procedure.
- Modular blades reduce overall costs by eliminating the need to replace the shaft.
- Calcar Wrench allows quick turnaround of instrument and intraoperative blade changes.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
1.5" Calcar Driver Complete Set (Hudson Drive End)	30-043-000
1.5" Calcar Driver Complete Set (Synthes Drive End)	30-045-000
1.5" Calcar Driver Complete Set (Tri-Shank Drive End)	30-041-000
1.8" Calcar Driver Complete Set (Hudson Drive End)	30-044-000
1.8" Calcar Driver Complete Set (Synthes Drive End)	30-046-000
1.8" Calcar Driver Complete Set (Tri-Shank Drive End)	30-042-000



**SYMMETRY Medical**<sup>®</sup> Inc.

486 West 350 North, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.7306  
www.symmetrymedical.com