


Charnley Retractor

The Charnley Retractor is designed for assistant-free retraction during minimally invasive hip and spine surgery.

BENEFITS

- Retracts on three sides (270 degrees).
- 33 different blade positions for maximum flexibility.
- Familiar attachment mechanism promotes quick, easy attachment and repositioning of blades.



Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Charnley Retractor, Complete Set	30-038-000
Standard Charnley Retractor Frame	13-045-008
Standard Charnley Retractor	13-045-001
Large Standard Charnley Retractor	13-045-002
Serrated Charnley Retractor Blade	13-045-003
Threaded Locking Screw	13-045-004



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com