

Patella Caliper

The Patella Caliper is useful for measuring the patella in total knee replacement. Restoring proper patellar thickness is essential for a successful outcome.

BENEFITS

- Easy-to-read calibrations range from 0 to 30mm.
- Ergonomic handles provide a secure grip.
- Easily cleaned

Contact your local Sales Manager for more information.



DESCRIPTION	PART #
Patellar Caliper	14-016-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com