

Screw Removal System

BENEFITS

- Using the locking thread design, the screw extractors remove screws by locking into the stripped hex head.
- Hex drivers are available with a cannulated shaft for easier removal of buried screws.
- Universal Extractor removes screws with heads partially or completely missing. The cone shaped head fully engages the remaining screw and optimizes the force needed for removal.
- Trephines are designed to fit over submerged screws for easy extraction with minimal bone loss.
- Universal Instrument Handle enables the surgeon to decide the appropriate rotational direction for a given application and wound exposure.
- The quick-connect release mechanism promotes quick interoperative exchange of instrument tips.



Symmetry Medical's Screw Removal System is designed to insert and remove solid and cannulated screws during trauma, hip, knee and arthroscopy procedures.

The intended uses for this system include removal of buried, stripped, and partial screws with broken screw heads.

The drive end (A/O) is designed for easy and quick engagement with the Universal Instrument Handle.

Contact your local Sales Manager for more information.

Screw Removal System

Product Descriptions

DESCRIPTION	PART #
Complete Screw Removal System	30-009-000
1.5mm Screw Extractor (A/O)	10-128-000
2.5mm Screw Extractor (A/O)	10-116-000
3.5mm Screw Extractor (A/O)	10-130-000
1.5mm Hex Driver (A/O)	10-117-000
2.5mm Hex Driver (A/O)	10-114-000
3.5mm Hex Driver (A/O)	10-115-000
4.0mm Hex Driver (A/O)	10-132-000
5.0mm Hex Driver (A/O)	10-133-000
Small Cruciform Screwdriver (A/O)	10-119-000
Mini Cruciform Screwdriver (A/O)	10-141-000
Straight Screwdriver (A/O)	10-120-000
Universal Extractor (A/O)	10-127-000
Small Extraction Bolt Body	10-127-003
Large Extraction Bolt Body	10-127-001
Fracture Reduction Pick	10-129-000
2.2mm I.D. Trephine (A/O)	10-121-000
3.2mm I.D. Trephine (A/O)	10-122-000
4.2mm I.D. Trephine (A/O)	10-123-000
4.7mm I.D. Trephine (A/O)	10-124-000
7.2mm I.D. Trephine (A/O)	10-125-000
Universal Driver Handle (A/O)	10-113-000
Radel® R Case	25-015-000



SYMMETRY Medical® Inc.

486 West 350 North, Warsaw, IN 46582, USA
 T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com