

# Patellar Tendon Staple

Symmetry Medical's Patellar Tendon Staple is designed for use during total knee replacement.

Securing the patellofemoral tendon during surgery can be cumbersome, while avoiding tibial tubercle avulsion. Normally, a K-Wire or a towel clip is used to hold the patellofemoral tendon in place.



## BENEFITS

- Offers an anatomic fit that accommodates all patient sizes.
- Wide band protects the tibial tubercle and patellar tendon.
- Angled staple legs increase pull-out resistance in poor quality bone.
- Designed to secure the patellofemoral tendon with two points of fixation during total knee replacement.
- Low profile design clears tibial resection guide and minimizes interference with saw blade.

DESCRIPTION	PART #
4 Degree Patellar Tendon Staple, Standard, Metal	14-025-000
8 Degree Patellar Tendon Staple, Large, Metal	14-024-000

Contact your local Sales Manager for more information.



**SYMMETRY Medical**<sup>®</sup> Inc.

486 West 350 North, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.7306  
www.symmetrymedical.com