

Universal Efficiency


Universal Instrument Handle

The Symmetry Medical Universal Instrument Handle is designed for use in two orientations: as a T-handle or a straight handle.

The single handle allows surgeons to decide which handle direction is most efficient and comfortable for a given application.

BENEFITS

- Reduce costs of entire set by eliminating handles attached to each instrument.
- Push button quick-connect release mechanism promotes quick intraoperative exchange of instrument tips.
- Ergonomic design reduces fatigue during lengthy procedures.



Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Universal Instrument Handle (Zimmer Drive End)	10-111-000
Universal Instrument Handle (Hudson Drive End)	10-112-000
Universal Instrument Handle (A/O Drive End)	10-113-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com