

Variable Angle Driver

The Variable Angle Acetabular Reamer Driver is designed to meet the various anatomies and patient sizes during Minimally Invasive Hip (MIH) Surgery.

The unique design incorporates a straight shaft with an adjustable head that can be set at three separate fixed points (per customer specifications) between 0 and 30 degrees.



BENEFITS

- Converts easily from *straight* driver to *angled* MIH driver for use during ALL Total Hip Replacement Surgeries.
- Three separate settings (0 to 30°) provide unmatched flexibility.
- Low-profile connection method provides improved visibility during MIH Surgery.
- Multiple driver/reamer options work with current reamer inventory.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Apex Connection (Zimmer/Tri-Shank End)	09-113-000
Bridgeback Connection (Zimmer/Tri-Shank End)	09-114-000
Apex Connection (Hudson End)	09-115-000
Bridgeback Connection (Hudson End)	09-116-000
Apex Connection (Synthes A/O End)	09-117-000
Bridgeback Connection (Synthes A/O End)	09-118-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com