

Minitainer


Protecting Your Investment

The Minitainer line of small to medium sized trays is specifically designed to help organize and protect delicate instruments during sterilization, transport and storage.

Minitainers are ideal for use with ophthalmic, ear/nose/throat, dental, cosmetic and microsurgery instruments.


Benefits

- Various sizes help customize procedural sets.
- Silicone mats gently secure instruments.
- Raised silicone (latex free) mats provide ventilation to facilitate drying.
- Clear Radel® R lids allow easy identification of instruments.



Description	Part #
2 x 7 Base, Lid and Mat, 7.12 x 1.92 x .73 inches (18 x 4.9 x 1.85 cm)	1-2700
3 x 6 Base, Lid and Mat, 5.84 x 2.5 x .7 inches (14.8 x 6.4 x 1.8 cm)	1-3600
3 x 7 Base, Lid and Mat, 7.1 x 2.86 x .73 inches (18 x 7.3 x 1.85 cm)	1-3700
4 x 8 Base, Lid and Mat, 7.5 x 4 x .75 inches (19 x 10.2 x 1.85 cm)	1-4800


**All dimensions are internal*



Minitainer Base, Lid and Silicone Bars

Description	Part #
2 x 7 Base and Lid with Silicone Bars for 1 Diamond Knife, 7.12 x 1.92 x .73 inches (18 x 4.9 x 1.85 cm)	1-2710
2 x 7 Base and Lid with Silicone Bars for 2 Diamond Knives, 7.12 x 1.92 x .73 inches (18 x 4.9 x 1.85 cm)	1-2720
3 x 7 Base and Lid with Silicone Bars for 3 Diamond Knives, 7.1 x 2.86 x .73 inches (18 x 7.3 x 1.85 cm)	1-3710
3 x 7 Base and Lid with 2 Silicone Bars for 8 Cannulae, 7.1 x 2.86 x .73 inches (18 x 7.3 x 1.85 cm)	1-3720


**All dimensions are internal*



Minitainer Base, Lid and Mat

Description	Part #
7 x 11 Single Level Base, Lid and Mat (Shallow), 10.5 x 6.25 x .62 inches (27 x 16 x 1.6 cm)	1-7100
7 x 11 Single Level Base, Lid and Mat (Deep), 10.5 x 6.25 x 1.36 inches (27 x 16 x 3.5 cm)	1-7110
7 x 11 Double Level Base, Lid and 2 Mats, 10.5 x 6.25 x 1.36 inches (27 x 16 x 3.5 cm)	1-7120


**All dimensions are internal*



Minitainer II Base, Lid and Mat

Description	Part #
6 x 8 Single Level Base, Lid and Mat, 5.65 x 7.25 x 1.2 inches (14.3 x 18.4 x 3 cm)	1-6800
7 x 11 Single Level Base, Lid and Mat, 6.18 x 10.25 x 1.2 inches (15.7 x 26.6 x 3 cm)	1-7200

**All dimensions are internal*



SYMMETRY Medical[®] Inc.

253 Abby Road, Manchester, NH 03103, USA
T: +1.603.647.7822 Fax: +1.603.647.7877

ZI de la Pilaterie, 14, Rue de la Ladrie, BP 165, 59653 Villeneuve d'Ascq Cedex, France
T: +33.3.28.09.94.54 Fax: +33.3.20.72.31.04