


Quick Launch System (QLS)



BENEFITS

- Over 25 silicone and Radel® brackets to accommodate a wide variety of instruments.
- Unique hole pattern allows brackets to be placed at different angles in 45 degree increments.
- All products are stocked for quick access and delivery.

Symmetry Medical offers QLS, which allows you to have what seems like a custom case configuration for just a small number of cases. These cases are ideal for prototype products and clinical trial instrument sets that may otherwise be damaged in handling and transport due to an inferior storage set.


You design the case configuration around your unique set of instruments using the three part case template. Symmetry PolyVac stocks all product cases, trays, lids and brackets for quick turnaround. This program allows safe and secure cases to be produced one at a time or in small quantities.

At Symmetry Medical, we are always looking for innovative ways to provide unique services to our customers. Contact your local Sales Manager for more information.


QLS Ordering Guide

STEP-BY-STEP GUIDE TO ORDERING YOUR QLS

Your Symmetry Sales Manager will be available to complete the template with you to assure accuracy and, ultimately, the best possible configuration for your components.



Step 1: Your Symmetry Representative will provide you with a three-part carbonless form that is the actual size of the case and tray.




Step 2: Arrange your products in a logical order on the template.



Step 3: Roughly trace around each part with a hard point pen. Press firmly.



Step 4: From a variety of brackets, select the brackets that best secure each component.



Step 5: Complete the information requested on the template including part numbers and quantity of each item required to complete your case.

Step 6: You can send the top copy of the template directly to Symmetry PolyVac. Keep one copy for yourself and make sure your Symmetry Sales Manager has one for reference.



253 Abby Road, Manchester, NH 03103, USA
T: +1.603.647.7822 Fax: +1.603.647.7877
www.symmetrymedical.com

ZI de la Pilaterie, 14, Rue de la Ladrie, BP 165, 59653 Villeneuve d'Ascq Cedex, France
T: +33.3.28.09.94.54 Fax: +33.3.20.72.31.04