



## Parallel Compressor & Distractor

The Symmetry Medical Parallel Compressor and Distractor are engineered to provide spinal compression and distraction between adjacent pedicle screws and hooks before final set-screw tightening.

With 30 degree angles, the narrow jaws improve visibility while fitting between closely placed pedicle screws. The parallel motion of the jaws maintain the instruments positions on the spinal rod, offering precise control.

Both Parallel Compressor and Parallel Distractor are designed to fit all spinal rod diameters. Instruments can be modified to meet customer requirements upon request.

- Parallel Compressor: 14-048-000
- Parallel Distractor: 14-049-000

Contact your local Sales Manager for more information.

# Engineered with PRECISION



486 West 350 North, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.7306  
[www.symmetrymedical.com](http://www.symmetrymedical.com)