



# HIP RETRACTORS

**Symmetry Medical's Hip Retractors** are designed for use in various surgical exposure to the hip, including minimally invasive approaches.

Typical Hohmann-type blade ends are combined with unique handle shapes for improved visualization and easy wound access. Each retractor thickness is optimized for ease of placement, visualization and stability.

- Typical Hohmann-type blade ends
- Retractor bend angles remove handles from wound site for improved visualization
- Optimized thickness helps wound integrity and stability
- Hip Retractor Complete System Part #: 30-034-000



**The trusted global source of innovative medical device solutions for today's needs and tomorrow's growth.**



**SYMMETRY Medical<sup>®</sup>** Inc.

3724 North State Road 15, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.4551  
[www.symmetrymedical.com](http://www.symmetrymedical.com)

## Hip Retractor Part Num

Hip Retractor, 25mm Blade	Part #: 13-023-000
Hip Retractor, Blunt Tip	Part #: 13-024-000
Offset Retractor, 32mm Blade	Part #: 13-025-000
Long Hohmann	Part #: 13-026-000
Long Hohmann, Curved, Short Blade	Part #: 13-027-000
Long Hohmann, Curved, Long Blade	Part #: 13-028-000
Long Hohmann, Narrow	Part #: 13-029-000
Small Hohmann	Part #: 13-030-000
Femoral Elevator	Part #: 13-031-000
Hybrid Case	Part #: 25-028-000

