



Minimally Invasive Knee (MIK) RETRACTOR SYSTEM

Our Minimally Invasive Knee (MIK) Retractor System is based on clinical use and surgeon feedback. The complete system features custom silicone handles for maximum function and comfort and unique shapes for MIK surgery.

The MIK Retractor System includes distinct retractors:

- Hohmann
- Superior
- Collateral
- Condylar
- Patellar
- Soft Tissue



The MIK Retractor System addresses several factors necessary for successful Total Knee Replacement (TKR) Surgery and Unicondylar Knee Replacement (UKR). Visibility during a minimally invasive approach is critical to ensure proper bone resections, ligament balancing and implant placement.

The MIK Retractors also protect soft tissue and ligaments. Several of the instruments are designed to provide retraction and also help prevent the saw blade from cutting through ligaments on the opposite side of the incision site.

Contact your local Sales Manager or visit www.symmetrymedical.com for more information.

The trusted global source of innovative medical device solutions for today's needs and tomorrow's growth.



SYMMETRY Medical[®]_{Inc.}

3724 North State Road 15, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.4551
www.symmetrymedical.com

Symmetry Medical's MIK Retractor System provides excellent visibility and unparalleled surgeon comfort during long TKR and UKR procedures.

MIK Retractor System Part Numbers and Descriptions



Small Hohmann Retractor

Part Number: 13-035-000

Used during initial dissection and designed to improve exposure of proximal tibia.



Large Hohmann Retractor

Part Number: 13-036-000

Wide Hohmann style blade for general retraction of large muscle groups.



Condylar Retractor

Part Number: 13-037-000

"J" curved blade used to retract soft tissue surrounding medial and lateral tibia during dissection.



Superior Retractor

Part Number: 13-038-000

Gentle radius of retractor designed to retract superior capsule without placing undue stress on skin incision.



Patellar Retractor

Part Number: 13-039-000

Windowed retractor blade captures patella, minimizing tension to patellar tendon and eliminating need for patellar eversion.



Collateral Retractor

Part Number: 13-040-000 (left), 13-041-000 (right)

Protect medial and lateral collateral ligaments from saw blade during proximal tibial resection.



Soft Tissue Retractor

Part Number: 13-042-000

Versatile soft tissue retractor used during exposure and initial bone resections.