



MODULAR CALCAR PLANER



The Symmetry Medical Modular Calcar Planer provides a safe, effective and accurate cutting instrument for the removal of hard calcar bone.

The patented design provides easy disassembly and blade change in a matter of seconds. The Calcar Wrench allows quick turnaround of instrument and intraoperative blade changes.

Surgeons will appreciate a constantly sharp blade producing repeatable, accurate cuts. The disposable blades provide sharp cuts for every procedure.

- Modular blades reduce overall costs by eliminating the need to replace the shaft
- The Calcar Wrench allows quick turnaround and intraoperative blade changes
- Disposable blades provide sharp cuts for every procedure

Contact your local Sales Manager or visit www.symmetrymedical.com for more information.



The trusted global source of innovative medical device solutions for today's needs and tomorrow's growth.



SYMMETRY Medical[®]

3724 North State Road 15, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.4551
www.symmetrymedical.com

Calcar Driver Part Numbers

* 1.5" Calcar Driver Complete Set Options

Complete Set with Hudson Drive End	Part #: 30-043-000
Complete Set with AO Drive End	Part #: 30-045-000
Complete Set with Zimmer/Tri-Shank Drive End	Part #: 30-041-000

* 1.8" Calcar Driver Complete Set Options

Complete Set with Hudson Drive End	Part #: 30-044-000
Complete Set with AO Drive End	Part #: 30-046-000
Complete Set with Zimmer/Tri-Shank Drive End	Part #: 30-042-000

Individual Instrument, Driver Only Options

1.5" Calcar Driver Only, Hudson Drive End	Part #: 11-103-000
1.5" Calcar Driver Only, AO Drive End	Part #: 11-104-000
1.5" Calcar Driver Only, Zimmer/Tri-Shank Drive End	Part #: 11-074-000
1.8" Calcar Driver Only, Hudson Drive End	Part #: 11-105-000
1.8" Calcar Driver Only, AO Drive End	Part #: 11-106-000
1.8" Calcar Driver Only, Zimmer/Tri-Shank Drive End	Part #: 11-076-000

Disposable Calcar Blades

1.5" Disposable Calcar Blade	Part #: 11-075-000
1.8" Disposable Calcar Blade	Part #: 11-077-000

Calcar Driver Wrench

Calcar Drive Wrench	Part #: 11-078-000
---------------------	--------------------



*Complete sets include: one modular calcar planer, one disposable calcar blade and one calcar driver wrench