



# Vertebral Preparation System

**Our Vertebral Preparation Anterior and Posterior Systems are** custom configured and ready to be rapidly deployed to your customers. Each system is complete with pre-defined selections for handles, shafts and tips. Expand your system with Symmetry Specialty Instruments for a customized solution.

The Vertebral Preparation Systems are built around our aesthetically and ergonomically designed silicone handles. Our handle options include one-handed or two-handed styles which can be customized in your choice of color mold with optional logo placement to provide company recognition.

- Pre-designed instruments help accelerate your product development cycle
- Complete Systems are pre-configured and ready to ship



**The trusted global source of innovative medical device solutions for today's needs and tomorrow's growth.**



**SYMMETRY Medical<sup>®</sup>** Inc.

3724 North State Road 15, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.4551  
[www.symmetrymedical.com](http://www.symmetrymedical.com)

## Vertebral Preparation Posterior System

### Complete 1-Handed Silicone Molded Instrument Handles System Part #: 30-056-000

Straight Osteotome 6mm	10-191-000
Angled Osteotome 6mm	10-192-000
Straight Cup Curette #4	10-193-000
Left Angled Cup Curette #4	10-194-000
Right Angled Cup Curette #4	10-195-000
Push Cup Curette #4	10-196-000
Pull Cup Curette #4	10-197-000
Straight Ring Curette 4mm	10-198-000
Angled Ring Curette 4mm	10-199-000
Bent Ring Curette 4mm	10-200-000
Rasp 7 x 16mm	10-201-000
Tamp 6mm	10-202-000
Straight Impactor	10-203-000
Angled Impactor	10-204-000
Mallet	10-205-000





## Vertebral Preparation Anterior System

**Complete 2-Handed Silicone Molded Instrument Handles  
System Part #: 30-055-000**

Cobb Elevator 25mm	10-206-000
Osteotome 12mm	10-207-000
Straight Cup Curette #3	10-208-000
Straight Cut Curette #6	10-209-000
Up-Angled Cup Curette #3	10-210-000
Up-Angled Cup Curette #6	10-211-000
Down-Angled Cup Curette #3	10-212-000
Down-Angled Cup Curette #6	10-213-000
Ring Curette 6mm	10-214-000
Ring Curette 9mm	10-215-000
Uterine Curette 12mm	10-216-000
Angled Rake Curette 12mm	10-217-000
Rasp 12mm	10-218-000
Tamp 9mm	10-219-000
Mallet	10-205-000



OPTIONS:

## Customize Your Vertebral Preparation System

Specify your tips, angles and shaft lengths


Select one or two-handed handle styles



Uterine Curettes



Penfield Elevators



Rake Curettes



Cup Curettes



Cobb Elevators



Round Graft Tamps



Rasps



Ring Curettes



Epstein Curettes

Contact your local Sales Manager or visit  
[www.symmetrymedical.com](http://www.symmetrymedical.com) for more information.



**SYMMETRY Medical<sup>®</sup>**

3724 North State Road 15, Warsaw, IN 46582, USA  
T: +1.574.267.8700 Fax: +1.574.267.4551  
[www.symmetrymedical.com](http://www.symmetrymedical.com)