

Bone Tamp


The Symmetry Medical Bone Tamp is designed specifically for use during ACL reconstruction.

The Bone Tamp securely compacts tunnel reamings into the remaining defect after harvest of the patellar plug for bone-tendon-bone autografts.

BENEFITS

- Ergonomic Radel® R handle design provides comfort and secure grip.
- Nine by 15mm tip fits typical defects.

Contact your local Sales Manager for more information.



DESCRIPTION	PART #
Bone Tamp with 9 x 15mm Tip, Ergonomic Radel R Handle	10-106-000



SYMMETRY Medical® Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com