

# Hybrid Reamers

Symmetry Medical's Hybrid Reamers are designed for use in minimally invasive or open total hip arthroplasty.

The unique three point connection method provides a secure connection between the driver and reamer while minimizing the incision needed for the procedure.

Hybrid Reamers are available as a complete set including drivers and cases.

## BENEFITS

- Apex connection improves vision into the acetabulum by eliminating connection "bars" on the back of the reamer.
- Incision length can be minimized as a result of reduced clearance needed for reamer/driver connection.
- Elongated tooth combined with Symmetry Medical's Smooth-Cut tooth design provides precise, full hemisphere reaming.
- Full size range in 1mm increments from 36mm to 80mm accommodates various patient anatomies.

Contact your local Sales Manager for more information.



# Hybrid Reamers

## Product Descriptions

DESCRIPTION	
36mm	16-021-036
37mm	16-021-037
38mm	16-021-038
39mm	16-021-039
40mm	16-021-040
41mm	16-021-041
42mm	16-021-042
43mm	16-021-043
44mm	16-021-044
45mm	16-021-045
46mm	16-021-046
47mm	16-021-047
48mm	16-021-048
49mm	16-021-049
50mm	16-021-050
51mm	16-021-051
52mm	16-021-052
53mm	16-021-053
54mm	16-021-054
55mm	16-021-055
56mm	16-021-056

DESCRIPTION	
57mm	16-021-057
58mm	16-021-058
59mm	16-021-059
60mm	16-021-060
61mm	16-021-061
62mm	16-021-062
63mm	16-021-063
64mm	16-021-064
65mm	16-021-065
66mm	16-021-066
67mm	16-021-067
68mm	16-021-068
69mm	16-021-069
70mm	16-021-070
71mm	16-021-071
72mm	16-021-072
73mm	16-021-073
74mm	16-021-074
75mm	16-021-075
76mm	16-021-076
77mm	16-021-077
78mm	16-021-078
79mm	16-021-079
80mm	16-021-080

