



Keeping a **LOW-PROFILE**



Parabolic Acetabular Reamers

The Symmetry Medical Parabolic Acetabular Reamer is the newest solution for reaming during Minimally Invasive Hip (MIH) Surgery.

The low-profile reamer inserts into the MIH incision smoothly, which helps minimize damage to the surrounding skin and soft tissue.

Reamer tooth geometry provides a smooth hemispherical acetabular surface as demonstrated by our standard Bridgeback Reamer while meeting the low-profile requirements of MIH exposure. The distinct reamer back design is intended to help reduce potential interference with retractors within the wound site.

Parabolic Reamers are compatible with all drivers featuring the patented Symmetry Bridgeback interface.

Contact your local Sales Manager or visit www.symmetrymedical.com for more information.



SYMMETRY Medical^{Inc.}

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com