

Modern minimally invasive surgical techniques for total hip replacement leave more muscle groups intact, resulting in difficulties when reducing the implants. The Namba Hip Slide from Symmetry Medical can be used as a lever to safely reduce a femoral implant into the acetabulum during total hip replacement.

While the Namba Hip Slide assists in reducing the articulation, the high-performance polymer of which it's made protects the implant surface – essential for ceramic heads. The Namba Hip Slide is radiolucent and can be used for femoral head diameters ranging from 22mm - 60mm. No special sterilization methods are required as the Namba Hip Slide is compatible with both steam and gas sterilization methods.

- 22 - 38mm heads Part #: 15-231-000
- 40 - 48mm heads Part #: 15-232-000
- 50 - 60mm Heads Part #: 15-233-000

