

Namba Hip Slide

Facilitates MIS hip replacement procedures



The new Namba Hip Slide gently slides the femoral head into the cup without scratching, essential for ceramic heads.

- Assists in reducing a femoral head trial and implant into the acetabulum during total hip surgery.
- Manufactured with high performance plastic material
- Radiolucent
- Compatible with steam or gas sterilization methods.

Namba Hip Slides accommodate 22 to 60mm diameter heads.

Contact your local Sales Manager or visit www.symmetrymedical.com for more information.

Namba Hip Slide (overall length 12"):

22 - 38mm heads: part # 15-231-000

40 - 48mm heads: part # 15-232-000

50 - 60mm heads: part # 15-233-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com

Cat. No. 5102-00 1M0209
Printed in USA. ©2009 Symmetry Medical Inc. All Rights Reserved.