


Symmetry Medical

Patellar Preparation System

The Patellar Preparation System accommodates all aspects of patellar replacement, regardless of technique or cementing process used.

This lightweight instrument reduces the amount of time dedicated to patella preparation during a total knee replacement.



BENEFITS

- Unique design promotes easy use.
- One-step *milling and drilling* process for single peg patella implant.
- Drop-in Drill Guide accommodates three peg patella implant.
- Single instrument serves as patellar reamer and cement clamp.
- System accommodates resurfacing and recessing techniques.
- Patella guides and graters are designed to accommodate specific sizes.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Patellar Reamer Clamp with Drive Shaft	14-026-000
Cement Clamp	14-026-059
Patellar Clamp Guide & Bushings	<i>Customer Specific</i>
Patellar Reamer Graters	<i>Customer Specific</i>
Without Pilot Hole	14-022-000-xx*
With Pilot Hole	14-026-061-xx*
Tri-Peg Drill Guide	<i>Customer Specific</i>
Single Peg Drill Guide	<i>Customer Specific</i>
Radel® R Case	25-025-000

*xx = reamer specific size



SYMMETRY Medical® Inc.

486 West 350 North, Warsaw, IN 46582, USA
 T: +1.574.267.8700 Fax: +1.574.267.7306
 www.symmetrymedical.com