

Symmetry Medical

Universal Hand Held Cases

Universal Hand Held Instrument Cases can accommodate a wide range of Hand Held instruments for use in arthroscopy, spinal, neurosurgery and microsurgery.

The Universal Hand Held Instrument Cases present the instruments in an upright position for easy access and identification on the sterile field.



BENEFITS

- Insert Carousel is removable to allow placement of additional instruments in the base (5-1150)
- Utility compartment and lid accommodate additional instruments (5-1110).
- Economical, durable and functional.

Contact your local Sales Manager for more information.

DESCRIPTION	PART #
Accommodates 8 Instruments Ranging from Rongeurs to Microsurgical Scissors with a Utility Space for Additional Instruments	5-1110
Accommodates up to 16 Instruments with Diameters up to 4mm, Additional Instruments can be Placed in Full Size Base	5-1150
Instrument Stand Accommodates up to 8 Small Hand Held Instruments in an Upright Position	5-1160



SYMMETRY Medical Inc.

253 Abby Road, Manchester, NH 03103, USA
T: +1.603.647.7822 Fax: +1.603.647.7877
www.symmetrymedical.com

ZI de la Pilaterie, 14, Rue de la Ladrie, BP 165, 59653 Villeneuve d'Ascq Cedex, France
T: +33.3.28.09.94.54 Fax: +33.3.20.72.31.04