

# 420 SS Acetabular Reamers

## BENEFITS

- Increased durability and longevity.
- Consistent cutting performance.
- Efficient preparation of the acetabulum.
- Exclusive material composition.
- Twenty percent more teeth than the competition for improved reaming of the acetabulum.
- Independent test data confirms 40 percent more durable than competitive reamers.<sup>1</sup>



Symmetry Medical reinforces its reputation for producing acetabular reamers with excellent cutting ability and longevity with the new 420 SS Acetabular Reamer.

Advances in material and teeth design guided the development of the 420 SS Acetabular Reamer with unmatched sharpness retention, durability and cutting consistency.<sup>2</sup>

Contact your local Sales Manager for more information.

<sup>1</sup> Sherry Lab, SMI 420SS Blast Reamer, July 25, 2005

<sup>2</sup> Dr. Susan Capps, PhD, Performance Evaluation of Symmetry 420SS Blast Acetabular Reamers, February 2006

# 420 SS Acetabular Reamers

## Product Descriptions

DESCRIPTION	
36mm	16-029-036
37mm	16-029-037
38mm	16-029-038
39mm	16-029-039
40mm	16-029-040
41mm	16-029-041
42mm	16-029-042
43mm	16-029-043
44mm	16-029-044
45mm	16-029-045
46mm	16-029-046
47mm	16-029-047
48mm	16-029-048
49mm	16-029-049
50mm	16-029-050
51mm	16-029-051
52mm	16-029-052
53mm	16-029-053
54mm	16-004-054
55mm	16-029-055
56mm	16-029-056

DESCRIPTION	
57mm	16-029-057
58mm	16-029-058
59mm	16-029-059
60mm	16-029-060
61mm	16-029-061
62mm	16-029-062
63mm	16-029-063
64mm	16-029-064
65mm	16-029-065
66mm	16-029-066
67mm	16-029-067
68mm	16-029-068
69mm	16-029-069
70mm	16-029-070
71mm	16-029-071
72mm	16-029-072
73mm	16-029-073
74mm	16-029-074
75mm	16-029-075
76mm	16-029-076
77mm	16-029-077
78mm	16-029-078
79mm	16-029-079
80mm	16-029-080



**SYMMETRY Medical**<sup>®</sup> Inc.

486 West 350 North, Warsaw, IN 46582, USA  
 T: +1.574.267.8700 Fax: +1.574.267.7306  
[www.symmetrymedical.com](http://www.symmetrymedical.com)